

RolaSlide®

Patient Transfer System

The Quintal RolaSlide® is a light weight patient transfer system that comprises a rolling slide board, 2 disposable slide sheets with handles and a pack of 10 disposable rollboard sleeves.

Patient transfers between two surfaces can now be undertaken easier, safer and more comfortably than with other lateral transfer systems. The system works perfectly in transferring patients from beds, x-ray tables, operating tables, shower trolleys and ambulance stretchers.


It reduces back injuries and prevents staff absence that can result from poorly implemented lateral patient transfer movements

The Disposable Rollboard sleeves keep the board free from contamination in use and reduce the need for frequent cleaning.



RolaSlide lateral patient transfer system

Exaggerated view



Ultra Low Friction Movement
Lightweight & Foldable

Fire Retardant
Anti-Static
Durable

RolaSlide's outstanding features ensure they are ideally suited for use in Hospitals, Care Homes and by the Emergency Services.

Call **0845 257 1717**
for further information

RolaSlide®

Patient Transfer System

PRODUCT DESCRIPTION

The Quintal RolaSlide is a lightweight, foldable, patient transfer system which transfers patients between a bed and x-ray table, operating table, shower trolley or ambulance stretcher. It efficiently replaces a patslide, draw sheet with one integrated system. the RolaGlide requires minimal effort to transfer a patient from one surface to another, and requires much less effort than a single slide sheet and patslide. The RolaSlide system is supplied with the Quintal roll board, 2 disposable slide sheets with handles and a pack of 10 disposable blue plastic rollboard sleeves. The roll board can be used on its own but then the patient will then have to be gently pushed across the surfaces by the carer.

PRODUCT SPECIFICATIONS

The Quintal RolaSlide comprises:

- 1 x RollBoard Patient Transfer Slide - 55cm wide x 177cm long (unfolded) and 55cm x 88cm (folded) &
 - 2 x NWF085190H disposable slidesheets with handles - 85cm wide x 190cm long
 - 1 x RSS01 pack of 10 Disposable RollBoard sleeves.
- Total WEIGHT 3.5 KG

SAFETY INFORMATION

- Before transferring patient, be sure both bed, stretcher, trolley, tables etc are securely locked
- Be sure both bed and stretcher, trolley tables are levelled when placing RolaSlide on them.
The maximum safety gap is 15cm between two surfaces for transfer..
- Do not sterilize the RolaSlide in high temperature
- When rolling the RolaSlide, be sure patients hair is positioned on the outer layer and away from the rolling to edge to prevent hair rolling up during operation
- Clean the RolaSlide using a surface disinfectant. Do not high temp sterile.
- The Safe Working Load is 150kg

Sizes shown are Nominal Sizes

RolaSlide Pricing

Product	Product Code	Comprises	Quantity in Pack	Price	
RolaSlide System	QRS01	RolaSlide Board & 2 x ProGlide slidesheets & 1xPkof 10 Rollboard sleeves	1 x Rollboard & 2 x slidesheets & 10 plastic sleeves	£232.00/sys	Equivalent to
Additional ProGlides	NWF085190H	Pack of 25 ProGlide Disposable Flatsheets With Handles	25/Pack	£43.00/pack	£1.72ea
Additional Sleeves	QRS02-DSL	Box of 50 Disposable Rollboard Sleeves	50/Box	£34.50/box	£0.69ea

**FREE DELIVERY
ON ALL ORDERS
OVER £500**

Quintal
HEALTHCARE

Unit 427A Birch Park
Thorp Arch Estate
Wetherby West Yorkshire LS23 7FG
Tel: +44 (0) 845 257 1717
Fax: +44 (0) 845 257 1818
www.quintalhealthcare.co.uk