


QProLateral™ Transfer Board

The **Qintal ProLateral Transfer Board** has been specifically designed to enhance the versatility of Lateral Transfers. Its unique design enables the board to either perfectly match the contours of a profiling bed or be used completely flat in a normal lateral transfer. In either configuration the board is designed to span a gap no greater than 150mm and this is clearly marked on the board

The ProLateral is particularly useful where a patient is either unable or afraid to lie completely flat during a lateral transfer.

Lightweight, with six ergonomically positioned handles the ProLateral board is easily placed and removed in any transfer situation. The handholds also ensure that the board will take up the minimum of space during storage.



Size 1500mm Long x 740mm Wide x 4.5 mm Thick

The ProLateral is lightweight, weighing less than 3.5kg and can be used either flat as a Patslide™ or curved with profiling bed transfers. The Anti-Slip strips on the reverse of the board provide maximum grip during transfer and it is easily wiped clean using soap & water, disinfectants or alcohol based solutions (up to 75%)

Also available as a rigid board- non folding same shape as QProLateral Designated QLat-slide

ProLateral's unique features ensure they are ideally suited for use in Care Homes, Hospitals and by the Emergency Services.

www.quintalhealthcare.co.uk

QProLateral

AET Transfer Board

The **Quintal QProLateral-AET** transfer board has been specifically designed to further extend the capabilities of the standard QProLateral Transfer Board. It has been designed narrower than the standard board in order to be easier to use within the restricted areas experienced by the Accident & Emergency Services and Operating Theatres.

With larger curved ends the AET board is also easier to position under the patient but its unique "hinged design" is still capable of perfectly matching the contours of a profiling stretcher or bed. Like the standard QProLateral the AET board, of course, can still be used completely flat in a normal lateral transfer like a "Patslide".



Size 1500mm Long x 590mm Wide x 4.5 mm Thick & Less than 3.4kg

In either configuration the board is designed to span a gap no greater than 150mm and this is clearly marked on the board

Anti-Slip strips on the reverse of the board provide maximum grip during transfer

The AET version of the QProLateral is invaluable when a patient, with breathing difficulties, is either unable or afraid to lie completely flat during a lateral transfer. It is lightweight, with six ergonomically positioned handles. It can be folded in half and the handholds also ensure that the board will take up the minimum of space during storage.



Also available as a rigid board-non folding same shape as QProLateral-AET
Designated
QLatslide-AET

Quintal
HEALTHCARE

Unit 427A Thorp Arch Estate
Wetherby West Yorkshire
LS23 7FG
Tel: +44 (0) 845 257 1717
Fax: +44 (0) 845 257 1818
www.quintalhealthcare.co.uk