


# Quintal ProLateral AET Board

The **Quintal ProLateral-AET Transfer Board** has been specifically designed to further extend the capabilities of the standard QProLateral Transfer Board. It has been designed narrower than the standard board in order to be easier to use within the restricted areas experienced by the Accident & Emergency Services and Operating Theatres.

With larger curved ends the AET board is also easier to position under the patient but its unique “hinged design” is still capable of perfectly matching the contours of a profiling stretcher or bed. Like the standard QProLateral the AET board, of course, can still be used completely flat in a normal lateral transfer like a “Patslide”.

In either configuration the board is designed to span a gap no greater than 150mm and this is clearly marked on the board



The AET version of the QProLateral is invaluable when a patient, with breathing difficulties, is either unable or afraid to lie completely flat during a lateral transfer. It is lightweight, with six ergonomically positioned handles. The handholds also ensure that the board will take up the minimum of space during storage.



- ◆ Lightweight less than 3.4kg
- ◆ Used either completely flat or curved to match a profiled trolley or bed
- ◆ 6 ergonomically placed hand holds make it easy to position and transport
- ◆ Can be folded in half for easy storage
- ◆ Clearly marked to show the board is designed to bridge a maximum gap of 150mm
- ◆ Anti-Slip strips on the reverse of the board provide maximum grip during transfer
- ◆ Easily wiped clean using soap & water, disinfectants or alcohol based solutions (up to 75%)

## Product Code

**QProLateral-AET** - Accident, Emergency & Theatre ProLateral Transfer Board  
Size - 1500mm x 590mm wide x 4.5 mm thick