


The **Quintal Range of QProLateral Transfer Boards** have been designed specifically to enhance the versatility of Lateral Transfers. Their unique design enables the board to either perfectly match the contours of a profiling bed or be used completely flat in a normal lateral transfer. The boards are printed with their optimum designed gap span of 150mm but with today's wider bedrails the boards can equally be used effectively with a gap up to a maximum of 180mm.

Lightweight, with six ergonomically positioned handles the QProLateral boards are easily positioned and removed in any transfer situation. The handholds also ensure that the boards take up the minimum of space during storage.

The curvature that the QProLateral boards can provide is particularly useful in circumstances where a patient may have breathing difficulties and is either unable or afraid to lie completely flat during a lateral transfer.



QProLateral Size
1500mm Long x 740mm Wide x 4.5 mm Thick



QProLat-AET Size
1500mm Long x 590mm Wide x 4.5 mm Thick

The QProLateral is designed to bridge the gap between two hospital beds. It is wider than a Patslide and provides more of a transfer surface with today's wider hospital beds. It is lightweight, weighing less than 4.5kg and is used either flat like a PatslideTM or curved to allow transfers between profiling beds or trolleys.

Anti-Slip strips on the reverse of the board ensure maximum grip during transfer. The QProLateral is easily wiped clean using soap & water, disinfectants or alcohol based solutions (up to 75%)

The board can be folded in half to ensure that it takes up the minimum of space during storage.

The QProLat-AET has been specifically designed to be narrower than the standard board. Within the restricted areas of A&E and Operating Theatres the larger curved ends of the AET board makes it easier to handle and position under the patient. It is lightweight, weighing less than 3.5kg and can also be used either flat like a PatslideTM or curved to allow transfers between profiling trolleys and theatre tables.


Anti-Slip strips on the reverse of the board provide maximum grip during transfer. The QProLat-AET is also easily wiped clean using soap & water, & propriety cleaning disinfectants or alcohol based solutions (to 75%)

The board can be folded in half to ensure that it takes up the minimum of space during storage.

ProLateral's unique features ensure they are ideally suited for use in Care Homes, Hospitals and by the Emergency Services.

The QProLateral Transfer Boards are also available in a non folding, completely flat rigid board format. These boards have the same shape as the curved boards but without the “hinges” that allow them to take up the curvature of profiling beds/trolleys. The boards are also printed with their optimum designed gap span of 150mm but with today’s wider bedrails the boards can equally be used effectively with a gap up to a maximum of 180mm.

Lightweight, with six ergonomically positioned handles the QProLateral boards are direct replacements to the Patslide. The handholds can be used to hang up the boards to take up the minimum of space during storage.



Qlatslide Size
1500mm Long x 740mm Wide x 4.5 mm Thick



Qlatslide-AET Size
1500mm Long x 590mm Wide x 4.5 mm Thick

QProLateral Range - Sizes & Pricing

Product	Product Code	Size (W x L) (Cm)	Type	Weight	Colour	Price
QProLateral	QProLateral	74 x 150	Hinged	4.5kg	White	£138.00
QProLat-AET1	QProLat-AET	59 x 150	Hinged	3.5kg	White	£119.00
Qlatslide	Qlatslide	74 x 150	Rigid	4.5kg	White	£138.00
Qlatslide-AET	Qlatslide-AET	59 x 150	Rigid	3.5kg	White	£119.00

FREE DELIVERY ON ALL ORDERS OVER £500

Sizes shown are Nominal Sizes